



WEIGAO MEDICAL INTERNATIONAL CO., LTD.

ADD: No. I Weigao Road, Hi-Tech Zone, Weihai, Shandong, China TEL: +86 631 5660699 FAX: +86 631 5622563 P.C.: 264210 E-mail: info@weigaogroup.com exp@weigaogroup.com Http://www.wegomedical.com/

WEGO


LEADING BRAND OF MEDICAL **DEVICES IN CHINA**











WG06-003-012 Drilled Plaster

Material: Drilled zine oxide with hot melt adhesive.


Waterproof or non-water proof is available

With drug or without is available

Color: White, Flesh.

Weight of adhesive: 60-65g/m2

Width: 2.5cm, 5cm, 7cm, 10cm, 12cm etc. Length: 1m, 5m etc.



WG06-003-013 Healing Capsicum Plaster

Capsicum plaster uses the full cotton cloth as the raw material, which was cocated on capsicum consisted zinc oxide adhesive. Capsicum plaster can make the partial skin generate physical reaction to accelerate blood fluxion. This plaster has drilled type and not-drilled type. All the products adopt the anti-sticky paper and film as the interleaver. Can be used by one piece or cut size conveniently and easily.



WG06-003-014 Sports Tape

Ideal for sport protection, water proof, permeability. different colors available.
Size available: 2.5cm×10y, 3.8cm×10y, 5cm×10y, 2.5cm×15y, 3.8×15y, 5cm×15y



WG06-003-015 Alcohol Prep Pad

Contents: 70% Isoprotyl Alcohol Non-woven Size: 65X30mm, 65X58mm I40XI80mm or other size per customer's demand. Type: I Ply, 2Ply, 4Ply, etc. Packet Size: 50X50mm, 60X85mm, etc.



WG06-003-016 PU Adhesive Tape

Material: Polyurethane film
Available in line shape or zig shape
Color: transparent
Size: 5X1000cm, 10X1000cm, 15X1000cm, 20X1000cm etc.



WG06-006-001 Blood Glucose Meter and Test Strips

High specific GOD, identify glucose accurately, avoid interference from other saccharides Automatically read the code of test strips HCT calibration automatically Auto power on/off when put in or pull out the test strip

Anti-interference formula Silver carbon high precision electrode Measuring method: Electrochemistry method



WG06-006-003/WG06-006-004 /WG06-006-005 Replaceable refill reusable insulin pen

Application fields: various liquid medicine injections, such as insulin,

GLP-I, aspart, glargine, liraglutide, etc Function: reusable, replaceable refill

Advantages: good feel / exquisite workmanship

Connection: Threaded connection for easy operation.

Accuracy range: The injection process is stable and the accuracy error is controlled at 2% -4% (international standard 5%). Injection range: 0-60IU, IIU increments, maximum 60IU (0.6ml)



WG06-006-002

- Silicone lubricant coating on insulin pen for smooth injection
- Easy to use, simply twisting on and off
- CE, ISO I 3485 certificate







WG06-004-001

One-piece Ostomy Pouch(Open Ostomy Pouch)

- •Specially designed hydrocolloid with high adhesion and no skin sensitivities
- •Soft and comfortable polymer film made bag body with opening design, easy to discharge
- Multiple types to satisfy clinical needs

WG06-004-002

Two-piece Ostomy Pouch(Open Ostomy Pouch)

- •Soft and comfortable polymer film made bag body with locking mechanism
- •Opening design eases discharging
- •Be able to apply with skin barrier to satisfy clinical needs

WG06-004-003 Ostomy Skin Barrier

•Higher erosion resistance, thinner flange, and better adhesion



Product Name	Picture	Туре	Specification	Opening/Flange Size
Two-Piece Ostomy Pouch		Drainable bag	21161	65mm/60mm
		Urostomy bag	22161	65mm/60mm



- Sterile: EO gas, non-toxic, non-pyrogenic.
- Size: 29G-32G
- Normal length: 4mm-12mm
- Color: different colors
- Diabetic care disposable products, for quick and convenient insulin injection, single used with insulin pen.

